

Digital Angiography System

DAR-9400f/9500f

IHE Integration Statement

REVISION HISTORY

REVISION	DATE	COMMENTS
Original	2007/11	Document is created.
A	2014/4	Apply to DAR-9500 Apply to Radiation Exposure Monitoring

Contents

1. Introduction	1
2. Connectathon	1
3. IHE Integration Statement.....	2

This page is intentionally left blank.

1. Introduction

This manual describes the DAR-9400f/9500f conformity to IHE (Integrating the Healthcare Enterprise).

IHE shows the ideal healthcare environment in the hospital.

To achieve the implementation along with IHE, Shimadzu have a willing to provide an ideal operation environment for our customers.

2. Connectathon

There is the test event called "Connectathon", is held in order to verify systems compliant with the IHE. Shimadzu had attended the North America Connectathon in 2014, and passed the prescribed test.

Note) DAR-9400f/9500f is not tested in IHE-J.

Please refer to the following site, to see the Connectathon test results.

<http://connectathon-results.ihe.net/>

Please refer to the following site, to see the IHE information of Shimadzu.

<http://www.shimadzu.com/products/medical/index.html>

IHE Integration Statement			
Vendor	Product Name	Version	Date
Shimadzu Corp	DAR-9400f, 9500f	Ver4.2, Ver5.2	2/27/14

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Patient Information Reconciliation	Acquisition Modality	None
Radiation Exposure Monitoring	Acquisition Modality	None
Scheduled Workflow	Image Display	None
Scheduled Workflow	Acquisition Modality	None

Internet address for vendors IHE information : <http://www.med.shimadzu.co.jp/>

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	http://www.shimadzu.com/products/medical/5iqj1d00000kvii.html

Link to general information on IHE : www.ihe.net